

#### NORTHWESTERN UNIVERSITY

### Addendum I: More about FDA The actual grant application

Pattanasak Mongkolwat PhD Research Associate Professor of Radiology

> David S. Channin MD Associate Professor of Radiology Chief, Imaging Informatics

> > www.radiology.northwestern.edu

### FDA Issues and Ideas

- "Base workstation" has clinical mode (FDA 510K)
- Switches to research mode on demand
  - Clearly marked to user
  - Automatic de-identification
  - Some users can only be in research mode
    - Site determination / IRB
- DICOM WG23 plug-ins: 2 kinds
  - Approved (Vendor/FDA)
  - WIP (works in progress)
    - Only available in research mode

#### FDA Issues and Ideas: Process

- Current NA-MIC/vTK/iTK builds/dashboards, etc.
  NECESSARY but not SUFFICIENT
- Must extend robust process back to design
  UML? Sparx Systems Enterprise Architect?
- Goal: Introduce process to FDA stepwise
  - Versus at end of development
- Goal: Freeze base workstation and take changes through FDA within 90 days

## The Grant

• 1R01-124379-01

"A Translation Workstation for Cancer Research Imaging"

- Submitted: January 17, 2006
- Pending Review: July 2006 (Yikes!)
- Assigned:
  - NCI
  - NIBIB
  - Roadmap Intitiatives
  - NIGMS

### The Team

- Northwestern
  - Pat Mongkolwat PhD
  - Alex Kogan
  - 2 x SW Engineer TBD
  - 2 x Grad Students TBD
  - 1 x Undergrad TBD
- Isomics
  - Stephen Pieper PhD
  - Alex Yarmarkovich PhD

# **Two Specific Aims**

- 1. Enhance NA-MIC's flagship software Slicer with:
  - 1. Clinical functionality
  - 2. Standards based interoperability
- 2. Test the null hypothesis:

There is no difference in the time it takes a diagnostic radiologist to perform a set of tasks related to interpreting a lung cancer imaging research study using a state-of-the-art commercial workstation and the time it takes to perform the same tasks on the enhanced Slicer workstation.

# **Clinical Functionality**

- Authoring tool for DICOM S/R objects
- Basic medical image display and manipulation
  - All 17 DICOM image objects (!)
  - Mapping DICOM objects directly into Slicer model (Isomics)
- DICOM WG23 support
- Hanging Protocols
- Annotation and Measurement
- Normal and Structured Reports (dictation out of scope)
- Security and de-identification

# **IHE Interoperability**

- IHE Consistent Presentation of Images
- IHE Reporting Workflow
- IHE Key Image Note
- IHE Teaching File and Clinical Trial Export
- IHE Evidence Document (for SIAM)
- IHE Audit Trail and Node Authentication

# Can we improve the research process?

- Give accession number
  to radiologist
- Retrieve study from worklist
- Retrieve historical case
- View current study
- View historical study
- Measure mass in new study

- Measure mass in old study
- Dictate report
- Export images from PACS
- Load study into 3D workstation
- Perform volumetric calculation of mass
- Dictate research report (from template)