




# Image Guided Therapy in Slicer3

Planning for Image Guided  
Neurosurgery

Danielle Pace, B.CmpH





# Acknowledgements

---



## Surgical Planning Lab, Harvard Medical School

Nobuhiko Hata, Steve Pieper, Ron Kikinis, Haiying Liu, Sonja Pujol, Randy Gollub, Rodrigo de Luis Garcia, Carl-Fredrik Westin



National Alliance for  
Medical Image  
Computing



Neuroimage  
Analysis Center



National Center for  
Image-Guided  
Therapy



Brain Science  
Foundation



# Acknowledgements

---

Thanks to

- **Alexandra Golby**, Department of Neurosurgery, Brigham and Women's Hospital
- **James Miller**, Visualization and Computer Vision, GE Research
- **Kilian Pohl**, Surgical Planning Laboratory, Department of Radiology, Brigham and Women's Hospital and Harvard Medical School
- **Takashi Suzuki**, Faculty of Advanced Techno-Surgery, Institute of Advanced Biomedical Engineering and Science, Tokyo Women's Medical University
- **Stephen Whalen**, Department of Neurosurgery, Brigham and Women's Hospital


for their helpful comments and suggestions



# Learning objectives

Following this tutorial, you will be able to perform all of the steps involved in preoperative planning for image guided neurosurgery within Slicer3, including:

- image registration
- model making
- DTI tractography





# Material

---

This course requires the installation of the Slicer3 software and the training datasets accessible at the following locations:

- Slicer3 software and building instructions:  
<http://www.slicer.org/pages/Downloads>
- Patient dataset, SPL-PNL brain atlas, and three pre-computed MRML scenes  
<http://wiki.na-mic.org/Wiki/index.php/IGT:ToolKit/Neurosurgical-Planning>


*Disclaimer: It is the responsibility of the user of 3D Slicer to comply with both the terms of the license and with the applicable laws, regulations and rules.*



# Prerequisites

- Data Loading and Visualization in Slicer3:



[http://wiki.na-mic.org/Wiki/index.php/Slicer:Workshops:Slicer3\\_Training](http://wiki.na-mic.org/Wiki/index.php/Slicer:Workshops:Slicer3_Training)





# 3D Slicer

- Integrates algorithms and utilities for medical image computing research and Image Guided Therapy into a single framework
- Is both an end-user application and a platform for research
- The precompiled program and the source code are both freely downloadable





# Image Guided Therapy (IGT) in Slicer3

---

Slicer3 has extensive support for IGT, including:

- Visualization
- Registration
- Segmentation
- Model making
- Diffusion Tensor Imaging
- Quantification
- Filtering
- Interfacing to imaging devices, trackers and medical robots




**Focus of  
this tutorial**





# The goal of neurosurgical planning

---



- Prior to surgery:
  - Integrate image information from multiple sources, including anatomical MRI, functional MRI and diffusion tensor imaging
  - Highlight structures of interest
  - Determine the best surgical approach




# Clinical Case - brain tumour resection

---

- Imaging showed a large lesion in the left frontal region of the brain, predicted to be a glioma (brain tumour originating from glial cells)
- Preoperative imaging included 3D SPGR, T2-weighted and FLAIR MRI imaging, language and motor functional MRI (fMRI) imaging, and diffusion tensor imaging (DTI)
- fMRI showed speech areas close to the lesion
- Surgical procedure: left frontal craniotomy and tumour resection

*See the [clinical\\_background](#) file within the patient dataset for more information*







# Tutorial datasets

---

- Three clinical datasets from a single patient with a large tumour in the left frontal region
- SPL-PNL brain atlas (based on a single healthy patient)




**SPGR MRI  
(anatomical)**



**Language  
fMRI**




**Diffusion  
Weighted  
Images**




**SPL-PNL  
brain atlas**




# Overview




1. **Loading and visualizing anatomical MRI data**




2. Incorporating fMRI data using image registration and thresholding




3. Creating a 3D model of the tumour volume



4. Predicting the locations of brain structures using image registration and a brain atlas



5. Incorporating brain fiber tractography from diffusion weighted images



6. Annotating the preoperative plan and saving the scene



# Anatomical imaging

---

- **Utility of anatomical imaging in IGT:**
  - Visualize brain structures of interest, such as the lesion
  - Perform measurements on structures of interest
  - Provide context for the rest of the scene



# Anatomical imaging

---


- **Steps involved in this section:**

Load anatomical MRI from patient dataset



# Load the anatomical MRI data


Open the  
Volumes  
module





# Load the anatomical MRI data

Click on  
“Select  
Volume File”








# Load the anatomical MRI data

Click on the first image in the 3D\_SPGR folder of the patient dataset, then click “Open”






# Load the anatomical MRI data

Keep the image origin set to “From File”

Click on “Apply” to load the image volume






# Load the anatomical MRI data

Click on the “link” button

Click on the “visibility” button

Click on the “fit to window” button


Click on the “center 3D view” button





# Load the anatomical MRI data

Open the Data module






# Load the anatomical MRI data

Click on the 000001.IMA node


Change the name of the MRML node to “anatomical-MRI” and press Enter






# Load the anatomical MRI data

Observe that the name of the node is changed in the MRML tree







# Overview




1. Loading and visualizing anatomical MRI data




**2. Incorporating fMRI data using image registration and thresholding**




3. Creating a 3D model of the tumour volume



4. Predicting the locations of brain structures using image registration and a brain atlas



5. Incorporating brain fiber tractography from diffusion weighted images




6. Annotating the preoperative plan and saving the scene



# Functional MRI (fMRI)

- **fMRI:**
  - Measures the blood oxygen level in each part of the brain while the patient performs a task, such as a speech or motor task
  - Statistical techniques are used determine which brain regions are active during the task
  - This statistical pre-processing has already been done using SPM

[\(http://www.fil.ion.ucl.ac.uk/spm/\)](http://www.fil.ion.ucl.ac.uk/spm/)



# SPM





# Functional MRI (fMRI)

---

- **Utility of functional MRI in IGT:**
  - Damage to regions of the brain important for language or movement could result in problems with speech, reading or movement
  - Knowing where these regions are allows us to modify our surgical plan so that we avoid them (as much as possible)



# Functional MRI (fMRI)

---

- **Steps involved in this section:**

Load language fMRI  
from patient dataset


Register language fMRI  
with anatomical MRI

Threshold language  
fMRI to find regions  
with high activation



# Load the language fMRI data


Open the  
Volumes  
module





# Load the language fMRI data


Click on  
“Select  
Volume File”





# Load the language fMRI data

Click on the first .hdr image in the Analyze folder within the language\_FMRI folder of the patient dataset, then click “Open”






# Load the language fMRI data

Set the Volume Name to languagefMRI

Keep the image origin set to "From File"

Click on "Apply"






# Load the language fMRI data

Set the foreground volume to languagefMRI

Set the background volume to none






# Load the language fMRI data

Set the background volume to anatomicalMRI

Click on the “fit to window” button

Click on the “center 3D view” button








# Load the language fMRI data


Scale the opacities to see both the foreground and background layers





# Load the language fMRI data


**Note that the image volumes are not aligned**






# Image-to-Image registration

---



The anatomical MRI and functional MRI image volumes are **not aligned**



The anatomical MRI and functional MRI image volumes **overlap**


**Image registration** aligns two images together with the goal of making the corresponding anatomy overlap



# Image-to-image registration

---

- Steps involved in image registration:




- An affine transformation estimates the rotation, translation, scale and shear needed to align the **moving** image with the **fixed** image
- Slicer's affine registration module uses **mutual information** to estimate the similarity between two images



# Initial manual transformation

Open the  
Transforms  
module






# Initial manual transformation

Create a new linear transform

The new transform is initialized to the identity matrix






# Initial manual transformation

Open the Data module

Change the name of the transform to “anatomical\_languagefMRI\_registrationTransform”






# Initial manual transformation

Drag the languagefMRI node under the transform node

Now any changes to the transform node are applied to the language fMRI image volume








# Initial manual transformation

Open the Transforms module

Adjust the translation components to roughly align the two image volumes in the axial plane






# Initial manual transformation

This initial transformation matrix “pushes” the language fMRI image onto the anatomical MRI image - the registration method expects the inverse

Click on the “Invert” button


Note that the images are no longer aligned





# Automatic affine registration


Open the  
Affine  
Registration  
module





# Automatic affine registration

Under  
Parameter  
Set, select  
“Create new  
CommandLine  
Module”






# Automatic affine registration

- Initial transform = anatomical\_language fMRI\_registration Transform
- Output transform = anatomical\_language fMRI\_registration Transform
- Fixed Image = anatomicalMRI
- Moving Image = languagefMRI
- Output Volume = Create New Volume

Click “Apply”






# Automatic affine registration

Set the foreground to languagefMRI, and the background to anatomicalMRI

Click on the “fit to window” button

Note that the images are not aligned






# Automatic affine registration

Open the Transforms module

Click "Invert"


The transform now transforms the moving image (languagefMRI) to align with the fixed image (anatomical MRI)





# Automatic affine registration

Click on the transform node to see that the matrix has changed








# Automatic affine registration

Open the Data module

The new volume in the MRML scene is the resampled moving image (the language fMRI)






# Automatic affine registration

Set the foreground to anatomical MRI and the background to Affine registration Volume1

Click on the “fit to window” button






# Automatic affine registration

Scale the opacities to see that the images are aligned


The new image is the moving image (language fMRI) transformed and resampled to align with the fixed image (anatomical MRI)





# Automatic affine registration

Right-click on the Affine registration Volume1 node and select "Delete node"





# Threshold image intensity


Open the Volumes module

Set the foreground to languagefMRI, and the background to anatomicalMRI

Click on the “fit to window” button

Set the active volume to languagefMRI


Expand the Display pane





# Threshold image intensity

Set the Color to FullRainbow






# Threshold image intensity


Set the threshold to Manual and adjust to highlight only high-intensity regions

Note that there are speech activation regions close to the tumour







# Overview




1. Loading and visualizing anatomical MRI data




2. Incorporating fMRI data using image registration and thresholding




**3. Creating a 3D model of the tumour volume**



4. Predicting the locations of brain structures using image registration and a brain atlas



5. Incorporating brain fiber tractography from diffusion weighted images



6. Annotating the preoperative plan and saving the scene






# Model Making in Slicer3

---

- **Utility of model making in IGT:**
  - View the structure as a complete unit (instead of slice by slice)
  - Perform measurements, such as volume measurements, that can be difficult to perform on the image volume itself





# Model Making in Slicer3

---

- **Steps involved in this section:**

Place a seed point on the tumour in the anatomical MRI

Segment the tumour volume in the anatomical MRI

Create the tumour model using Slicer3's ModelMaker


Change the tumour model's appearance

Save the tumour model



# Place the seed point


Click on the Background layer button to show the anatomical MRI alone





# Place the seed point


Open the  
Fiducials  
module





# Place the seed point


Create a new Fiducial List





# Place the seed point


The fiducial will show up as a pink starburst





# Place the seed point

Set the mouse mode to “place items”






# Place the seed point

Slice through the image volume until you see the tumour in all three views

Click once to place the seed








# Place the seed point


The fiducial can also be seen in the 3D Viewer





# Place the seed point


Set the mouse mode to “transform view”





# Segment the tumour volume

Open the Simple Region Growing module






# Segment the tumour volume

- Parameter set = Create new CommandLine Module
- # of iterations = 1
- Seeds = FiducialList1
- Input Volume = anatomicalMRI
- Output Volume = Create New Volume

Click “Apply”

This step may take a few minutes to run - wait until the status says “Completed”






# Segment the tumour volume

Set the background to anatomicalMRI

Click on the “fit to window” button


The tumour is segmented by the label layer called “Simple Region growing Volume1” (shown in purple)





# Create the tumour model


Open the Model Maker module





# Create the tumour model

- Parameter set = Create new CommandLine Module
- Input Volume = Simple region growing Volume1
- Models = Create New ModelHierarchy
- Model Prefix = Tumour
- Labels = 255 (scroll over the purple tumour and note the “Lb” value)






# Create the tumour model

Click “Apply”

This step may take a few minutes to run - wait until the status says “Completed”








# Create the tumour model

Slice through the image volume to expose the tumour model in the 3D Viewer






# Create the tumour model

Open the Data module

Delete the FiducialList and Simple region growing Volume by right-clicking and selecting "Delete Node"






# Create the tumour model

Note that the model now appears in the MRML tree


Change the name of the model to "tumourModel"





# Change the model's appearance


Open the  
Models  
module





# Change the model's appearance

Select the tumour model






# Change the model's appearance

Set the model's opacity to 0.6


Click on the box next to "Set Color" to change the colour of the model to blue





# View the anatomical MRI, fMRI and tumour model


Click on the Foreground layer button to show the language fMRI data





# Save the tumour model

In the  
“Save”  
pane, click  
“Save  
Model”








# Save the tumour model


Name the tumour model  
tumour.vtk

Click  
“Save”







# Overview




1. Loading and visualizing anatomical MRI data




2. Incorporating fMRI data using image registration and thresholding




3. Creating a 3D model of the tumour volume



4. **Predicting the locations of brain structures using image registration and a brain atlas**



5. Incorporating brain fiber tractography from diffusion weighted images



6. Annotating the preoperative plan and saving the scene



# The SPL-PNL Brain Atlas

---

[http://www.na-mic.org/pages/Special:PubDB\\_View?  
dSPACEID=1265](http://www.na-mic.org/pages/Special:PubDB_View?dSPACEID=1265)



- Includes:
    - Anatomical MRI
    - Label maps
    - >160 models
- from a healthy volunteer



# Incorporating a brain atlas

---

- **Utility of an atlas in IGT:**
  - Estimate the positions, orientations and shapes of important brain structures without having to segment each one in the patient dataset



# Incorporating a brain atlas

---


- **Components of the brain atlas used:**
  - Anatomical MRI image
  - Models of segmented brain structures
- **Using a brain atlas to approximate where brain structures are in the patient:**
  - Register the atlas's MRI with the patient's MRI to determine the transformation between the atlas and the patient
  - Use the same transformation to transform the atlas models so that they overlap with the patient's MRI



# Incorporating a brain atlas

---

- **Steps involved in this section:**




- The precomputed MRML scene contains a copy of the patient's anatomical MRI, the atlas's anatomical MRI, and three models:
  - Left and right optic tracts
  - Left anterior superior temporal gyrus
- The anatomical MRIs of the patient and the atlas were registered with an initial manual transform followed by an automatic affine registration



# Import the patient-atlas registration transform scene

Select File -> Import Scene


Note: Make sure that you select “Import Scene” and not “Load Scene”, as “Load Scene” will delete your work so far





# Import the patient-atlas registration transform scene

Select  
“anatomical  
MRI  
Registered  
ToAtlas.mrml”,  
then click  
“Open”








# Inspect the MRML scene

Open the Data module

Set the foreground to atlasMRI and the background to anatomicalMRI

Click on the “fit to window” button






# Inspect the MRML scene

Inspect the atlas MRML scene - 6 new nodes

- atlasMRI = anatomical MRI from the atlas
- optic\_tract\_L.vtk = model of the left optic tract
- optic\_tract\_R.vtk = model of the right optic tract
- Superior\_temporal\_gyrus\_ant\_L.vtk = model of the left anterior superior temporal gyrus






# Inspect the MRML scene

Inspect the atlas MRML scene - 6 new nodes


- anatomical MRI = copy of the patient's anatomical MRI
- anatomicalMRI\_atlas\_registration Transformation = the transformation between the patient's anatomical MRI and the atlas's MRI





# Inspect the MRML scene


Delete the duplicate anatomicalMRI node by right-clicking and selecting “Delete Node”





# Inspect the patient-atlas registration

Scale the opacities to see how well the atlas MRI and patient anatomical MRI image volumes are aligned





# Inspect the patient-atlas registration

---

- **Important note:** Since an affine registration was used, the brain structure models from the atlas may not match the anatomy of the patient:
  - Affine registration does not account for differences in brain shape
  - The patient's tumour will deform brain structures compared to the healthy atlas subject
- Thus the brain structure models give only a rough idea of their positions in the patient
- Non-linear registration could be used to partially remedy this - ex. B-Spline registration in Slicer3




# Inspect the entire scene so far

Set the foreground to languagefMRI


Click on the Foreground layer button to show the language fMRI data

Slice through the image volume to get a good view







# Overview




1. Loading and visualizing anatomical MRI data




2. Incorporating fMRI data using image registration and thresholding




3. Creating a 3D model of the tumour volume



4. Predicting the locations of brain structures using image registration and a brain atlas



**5. Incorporating brain fiber tractography from diffusion weighted images**




6. Annotating the preoperative plan and saving the scene






# Diffusion Tensor Imaging (DTI)

## Diffusion Sensitizing Gradients



## Diffusion Weighted Images






- Each diffusion weighted image shows the diffusion of water in the direction of the diffusion sensitizing gradient



# Diffusion Tensor Imaging (DTI)

---

- Acquire diffusion weighted images (DWI) and baseline images
  - Estimate a tensor at each voxel to create a Diffusion Weighted image (DWI)
    - Calculate a scalar “diffusion anisotropy index”, such as fractional anisotropy, to quantify the diffusion’s anisotropy
      - Perform tractography (fiber tracking)





# Diffusion Tensor Imaging (DTI)


---

- **Utility of Diffusion Tensor Imaging in IGT:**
  - Major tracts that run between functionally important regions of the brain must remain intact to prevent side effects for the patient
  - Knowing where these tracts are allows us to modify our surgical plan so that we avoid them (as much as possible)



# Diffusion Tensor Imaging (DTI)

- **Steps involved in this section:**




Registration using a precomputed MRML scene

- The precomputed MRML scene contains a copy of the patient's anatomical MRI, a copy of the baseline image of the patient's DWI, and the transformation that registers them together



# Load the DWI data

Open the  
DICOM DWI  
Loader module






# Load the DWI data

Create a new  
CommandLine  
Module


Click on the  
folder icon to  
set the  
DICOM  
directory





# Load the DWI data

Click on the DTI folder within the patient dataset, then click “OK”






# Load the DWI data

Set the dwiimage to “Create New Diffusion Weighted Volume”

Click “Apply”

This step may take a few minutes. Wait until Status says “Completed”








# Load the DWI data

Open the Data module


Rename the new DICOM DWI Loader Volume node to “DTI” by selecting it and changing the name in the MRML Node Inspector pane





# Load the DWI data

Turn the visibility off for the four models in the scene (tumour, left and right optic tracts, and left anterior superior temporal gyrus) by right-clicking and selecting “Toggle Visibility”






# Inspect the DWI data

Open the Volumes module

Open the Display pane

Set the active volume to the DTI volume






# Inspect the DWI data

Note that the Display pane now allows you to select the DWI volume

Set the background to DTI

Click on the Background layer button

Click on the “fit to window” button






# Inspect the DWI data

Scroll through the 36 volumes in the set:

-0 to 4 are baseline images

-5 to 35 are diffusion weighted images






# Create the tensors

Note that the Diffusion Editor pane within the volumes module is activated when diffusion weighted images are loaded

Expand the Diffusion Editor pane and note the gradients used to acquire the images






# Create the tensors

Expand the “Testing (Tensor Estimation & Glyphs & Tractography Seeding)” pane

You would click “Run” to compute the tensors, but will upload a MRML scene of precomputed tensors because computing tensors takes a while






# Create the tensors

Open the Data module

Select File -> Import Scene

Note: Make sure that you select “Import Scene” and not “Load Scene”, as “Load Scene” will delete your work so far








# Create the tensors

Select  
“precomputed  
Tensors.mrml”,  
then click  
“Open”






# Create the tensors

3 new nodes:

- **Baseline Node:** image without diffusion weighting
- **Threshold Mask:** mask excluding the background
- **TensorNode** = the tensors

Set the background to the tensor node






# Create the tensors

Open the Volumes module

Set the active volume to the diffusion tensor node


The fractional anisotropy is displayed by default





# View the tensors

Set the Window/Level to Manual and adjust the levels to invert the image so that the tensors will be easier to see






# View the tensors

Scroll down and turn the glyphs on for all three views

Adjust the scale factor and the spacing


The line's direction is the direction of the tensor's main eigenvector, while the colour shows the fractional anisotropy





# View the tensors


Change the  
Glyph Type to  
“Ellipsoids”





# View the tensors

Change the  
Glyph Type  
back to “Lines”






# Place the tractography seed points

Open the  
Fiducials  
module

Create a new  
Fiducial List








# Place the tractography seed points


Set the glyph  
type to  
Sphere3D





# Place the tractography seed points

Set the mouse mode to “place items”






# Place the tractography seed points

Click to place approximately ten fiducial points along the corpus callosum


Set the mouse mode to “transform view”





# Perform tractography (fiber tracking)


Open the tractography “Fiducial Seeding” module





# Perform tractography (fiber tracking)

- DTI volume = Diffusion Editor\_1\_TensorNode
- Fiducial list = FiducialList2
- Output FiberBundle Node = Create New Fiber Bundle






# Perform tractography (fiber tracking)

Open the Volumes module

Turn the glyphs off for all three views

Click on the “visibility” button to turn off the slice visibility






# Perform tractography (fiber tracking)

One tract is generated for each fiducial.

The tract colour is the fractional anisotropy by default.






# Perform tractography (fiber tracking)

Open the Fiducials module

Turn off the visibility of the fiducials








# Import the anatomical-DTI registration transform

Select File -> Import Scene


Note: Make sure that you select “Import Scene” and not “Load Scene”, as “Load Scene” will delete your work so far





# Import the anatomical-DTI registration transform

Select  
“anatomical  
MRI  
Registered  
ToDTI.mrml”,  
then click  
“Open”






# Inspect the MRML scene

Open the Data module

Inspect the atlas MRML scene - 3 new nodes


- anatomicalMRI = copy of the patient's anatomical MRI
- BaselineNode = copy of the patient's DTI baseline node
- anatomicalMRI\_DTI\_registrationTransform = aligns the DTI images with the anatomical image





# Inspect the MRML scene


Delete the duplicate anatomicalMRI node by right-clicking and selecting “Delete Node”





# Inspect the MRML scene


Delete the duplicate BaselineNode node by right-clicking and selecting “Delete Node”





# Inspect the MRML scene

Drag the three diffusion editor nodes, the fiducial list used to seed the tractography and the fiber bundle node under the anatomicalMRI\_DTI\_registrationTransform






# Inspect the anatomical MRI-DTI registration

Set the foreground to the baseline DTI node and the background to anatomicalMRI

Scale the opacities to see how well the atlas MRI and patient anatomical MRI image volumes are aligned






# Inspect the anatomical MRI-DTI registration

Set the foreground to the tensor node and the background to anatomicalMRI

Scale the opacities to see how well the atlas MRI and patient anatomical MRI image volumes are aligned








# Inspect the entire scene so far

Set the foreground to languagefMRI

Click on the “visibility” button


Click on the Foreground layer button to show the language fMRI data





# Inspect the entire scene so far


Slice through the image volume to get a good view





# Inspect the entire scene so far


Turn the visibility on for the four models in the scene (tumour, left and right optic tracts, and left anterior superior temporal gyrus) by right-clicking and selecting “Toggle Visibility”






# Inspect the entire scene so far

The scene now contains anatomical MRI data, functional MRI data, a segmented tumour volume, estimates of other brain structures and brain fiber tractography







# Overview




1. Loading and visualizing anatomical MRI data




2. Incorporating fMRI data using image registration and thresholding




3. Creating a 3D model of the tumour volume



4. Predicting the locations of brain structures using image registration and a brain atlas



5. Incorporating brain fiber tractography from diffusion weighted images




**6. Annotating the preoperative plan and saving the scene**



# Plan annotation

---

- **Utility of plan annotation in IGT:**
  - Surgeons can mark important points on the plan, such as the planned surgical access point, so that their notes can be brought into the operating room





# Plan annotation

---

- **Steps involved in this section:**

Mark the entry point and the tumour access point on the plan

Measure the distance between the two points

Save the scene



# Annotate the plan

Open the  
Fiducials  
module

Create a new  
Fiducial List

The screenshot shows the 3D Slicer 3.0 interface. The Fiducials module is active, and a new fiducial list named 'FiducialList1' has been created. The 3D view shows a brain model with several red fiducial markers. The 2D views (Axial, Sagittal, Coronal) also show the brain with the markers. The 'Create New Fiducial List' dialog is open, and the 'Fiducial List Select' dropdown is set to 'FiducialList1'. The 'Text Scale' is set to 4.5, and the 'Symbol Scale' is set to 10. The 'Opacity' is set to 1. The 'Set Color' and 'Set Selected Color' buttons are visible. The 'Set Opacity Type' is set to 'Sphere3D'. The table below shows the list of fiducial markers:


Name	Selected	X	Y
FiducialList1-P1	<input checked="" type="checkbox"/>	0.00...	-30.4...
FiducialList1-P2	<input checked="" type="checkbox"/>	0.00...	-16.7...
FiducialList1-P3	<input checked="" type="checkbox"/>	0.00...	-4.08...
FiducialList1-P4	<input checked="" type="checkbox"/>	0.00...	4.99...
FiducialList1-P5	<input checked="" type="checkbox"/>	0.00...	14.9...
FiducialList1-P6	<input checked="" type="checkbox"/>	0.00...	24.0...
FiducialList1-P7	<input checked="" type="checkbox"/>	0.00...	31.3...





# Annotate the plan


Set the glyph  
type to  
Sphere3D





# Annotate the plan


Slice through the sagittal plane until it roughly bisects the tumour volume





# Annotate the plan


Adjust the 3D Viewer to get a view from the right





# Annotate the plan


Set the mouse mode to “place items”





# Annotate the plan


Place a fiducial on the sagittal plane representing the planned entry point for the left frontal craniotomy





# Annotate the plan


Double-click on the name of the fiducial and change it to "entryPoint"





# Annotate the plan


Place a fiducial on the sagittal plane representing the planned access point to the tumour





# Annotate the plan

Double-click on the name of the fiducial and change it to "tumour Access"








# Distance measurements

Scroll down, and note that the distance between the two fiducials is shown


Change the mouse mode to “transform view”





# Inspect the final scene


Slice through the image volume to get a good view





# Inspect the final scene


Set the layout to the “3D only layout”





# Inspect the final scene


The result of neurosurgical planning





# Save the scene

Select File ->  
Save






# Save the scene

Click on the “Scene File” icon to name the .mrml scene and to specify its directory

Click “Save”





# Overview

---

- **In this tutorial, you learned how to**
  - Register image volumes together
  - Segment and build models of structures of interest
  - Load diffusion weighted images, calculate tensors and perform fiber tractography

**in order to build a preoperative neurosurgical plan incorporating**

- anatomical MRI
- functional MRI
- a model of the tumour
- brain structure models from a registered atlas
- fiber tractography
- annotated fiducial points



# Conclusions

---

- Slicer3 has extensive support for Image Guided Therapy
- Slicer3 is free open-source software that allows IGT researchers to share algorithms and work within a common framework





# For more information...

---

## Image registration:

- Read about the mutual information similarity measure here:

Wells, S. *et al.*, Multi-modal volume registration by maximization of mutual information. *Medical Image Analysis*, 1(1):35-51 (1996).

Collignon, A. *et al.*, Automated multimodality image registration based on information theory. *Proceedings of the International Conference on Information Processing in Medical Imaging*, 263-274 (1995).



# For more information...

---

## Model Making:

- Slicer3 uses the marching-cubes algorithm to generate models - read about it here:

Lorensen, W.E. and Cline, H.E. Marching cubes: a high resolution 3D surface construction algorithm. ACM SIGGRAPH Computer Graphics 21(4):163-169, 1987.



# For more information...

---

## Diffusion Tensor Imaging:

- A review on DTI:

Le Bihan, D. *et al.*, Diffusion Tensor Imaging: concepts and applications. *Journal of Magnetic Resonance Imaging*, 13:534-546 (2001)

- A review on tractography:

Mori, S. and van Zijl, P.C.M., Fiber tracking: principles and strategies - a technical review. *NMR in Biomedicine*, 15:468-480 (2002)