



Diffusion MRI Analysis

Sonia Pujol, Ph.D.

Surgical Planning Laboratory,
Harvard Medical School

spujol@bwh.harvard.edu

Brain Anatomy




- White matter ~45% of the brain
- Myelinated nerve fibers (~ 10 μm axon diameter)

White Matter Exploration



Jules Joseph Dejerine (*Anatomie des centres nerveux* (Paris, 1890-1901): Atlas of Neuroanatomy based on myelin stained preparation


White Matter Exploration



First non-invasive
window on the
organization of brain
white matter
pathways *in-vivo*


Tutorial Outline

This tutorial is an introduction to the fundamentals of Diffusion MRI analysis, from the estimation of diffusion tensors to the interactive 3D visualization of fiber tracts in Slicer4.2



Tutorial dataset

The tutorial dataset DiffusionMRI_tutorialData is a Diffusion Weighted MR scan of the brain acquired with 41 gradient directions and one baseline.



The dataset is available on the Slicer Training Compendium (www.slicer.org)

Tutorial software




The tutorial uses the
3DSlicer version 4.2
software available at
www.slicer.org

Disclaimer

It is the responsibility of the user of 3DSlicer to comply with both the terms of the license and with the applicable laws, regulations and rules. Slicer is a tool for research, and is not FDA approved.

3DSlicer



3D Slicer is a multi-institution effort supported by the National Institutes of Health.


- An **end-user application** for image analysis
- An **open-source environment** for software development
- A software platform that is both **easy to use** for clinical researchers and **easy to extend** for programmers

Learning Objectives


Following this tutorial, you'll be able to

- 1) Estimate a tensor volume from a set of Diffusion Weighted Images
- 2) Understand the shape and size of the diffusion ellipsoid
- 3) Reconstruct DTI tracts from a pre-defined region of interest
- 4) Interactively visualize DTI tracts seeded from a fiducial


MR Diffusion Analysis Pipeline




DWI
Acquisition




Tensor
Calculation



Scalar
Maps




3D
Visualization




Part 1: From DWI images to Tensors

Understanding the DWI dataset




The Diffusion Weighted Imaging (DWI) dataset is composed of 1 volume acquired without diffusion-sensitizing gradient, and 41 volumes acquired with 41 different diffusion-sensitizing gradient directions.

Start the Slicer Software




None RAS: (125.0, -125.0, 1.0),

Start the Slicer software




Slicer displays the list of 103 modules in alphabetical order.


Loading the DWI dataset




Loading the DWI dataset




Loading the DWI dataset




Loading the DWI dataset



Loading the DWI dataset



Adjusting Window and Level



The screenshot displays the 3D Slicer software interface. The top menu bar shows 'Modules' and 'Volumes'. The left sidebar contains a list of modules, with 'Volumes' selected. The main window shows a 3D view of a brain MRI scan with a purple wireframe box indicating the current slice. A red arrow points to the 'Threshold' slider in the 'Volumes' module panel, which is currently set to 'Off'. The 'W/L' (Window/Level) settings are visible, showing 'W: 532' and 'L: 272'. The 'Histogram' button is also present. The bottom of the interface shows three orthogonal views (axial, sagittal, and coronal) of the brain scan, with color-coded axes (R, Y, G) and a scale bar.

Select the module Volumes from the modules menu.

Adjust the window and display of the baseline image using the W/L slider

Exploring the DWI dataset

The baseline image corresponds to the DWI Component #0.

Select the DWI component #10, which corresponds to the 10th diffusion sensitizing gradient

W: 1478 Manual W/L L: 529

-1000 4044


0.00 4040.00

Histogram


Data Probe

L
F
B


Exploring the DWI dataset



Exploring the DWI dataset



Exploring the DWI dataset




Exploring the DWI dataset

Left click on the pin button in the top left corner of the red viewer to display the slice menu.


Click on the 'links' icon to link all three viewers, and click on the 'fit image to window icon'.

L
F
B

Exploring the DWI dataset




Diffusion Tensor Estimation




Select the module **DWI to DTI Estimation** in the modules menu:

- select the Input DWI volume 'dwi'
- select Output DTI Volume 'Create New Diffusion Tensor Volume', and rename it 'dti'
- select Output Baseline Volume 'Create new Volume', and rename it 'baseline'
- select the Estimation Method 'WLS' (Weighted Least Squares) and click on Apply.


Diffusion Tensor Estimation



Diffusion Tensor Estimation



Diffusion Tensor Data



$$S_i = S_0 e^{-b \hat{g}_i^T \underline{D} \hat{g}_i}$$


Stejskal-Tanner equation (1965)

$$\underline{\mathbf{D}} = \begin{bmatrix} D_{xx} & D_{xy} & D_{xz} \\ D_{yx} & D_{yy} & D_{yz} \\ D_{zx} & D_{zy} & D_{zz} \end{bmatrix}$$


The diffusion tensor $\underline{\mathbf{D}}$ in the voxel (I,J,K) is a 3x3 symmetric matrix.

Diffusion Tensor

- The diffusion tensor \underline{D} in the voxel (I,J,K) can be visualized as an ellipsoid, with the eigenvectors indicating the directions of the principal axes, and the square root of the eigenvalues defining the ellipsoidal radii.
- Scalar maps can be derived from the rotationally invariant eigenvalues λ_1 , λ_2 , λ_3 to characterize the size and shape of the diffusion tensor.




Diffusion Tensor Shape




$$\lambda_1 = \lambda_2 = \lambda_3$$

Isotropic media
(CSF, gray matter)




$$\lambda_1 \gg \lambda_2, \lambda_3$$

Anisotropic media
(white matter)



$$\lambda_1 \sim \lambda_2 \gg \lambda_3$$

Exploring the Diffusion Tensor Data



Corpus Callosum

The corpus callosum is a broad thick bundle of dense myelinated fibers that connect the left and right hemisphere. It is the largest white matter structure in the brain





Image from Gray's Anatomy

Exploring the Diffusion Tensor Data




Characterizing the Size of the tensor: Trace


$$\text{Trace}(D) = \lambda_1 + \lambda_2 + \lambda_3$$

- $\text{Trace}(D)$ is intrinsic to the tissue and is independent of fiber orientation, and diffusion sensitizing gradient directions
- $\text{Trace}(D)$ is a clinically relevant parameter for monitoring stroke and neurological condition (degree of structural coherence in tissue)
- $\text{Trace}(D)$ is useful to characterize the size of the diffusion ellipsoid


Characterizing the Size of the tensor: Trace




Trace




Trace



Trace




Move the mouse cursor in the 2D view, and observe the values of the trace in the corpus callosum and in the adjacent gray matter.






Trace

Note how the Trace values are fairly uniform in both white and gray matter, even if the tissues are different in structure.




Scalar Maps: Fractional Anisotropy

$$FA(D) = \frac{\sqrt{(\lambda_1 - \lambda_2)^2 + (\lambda_1 - \lambda_3)^2 + (\lambda_2 - \lambda_3)^2}}{\sqrt{2} \sqrt{\lambda_1^2 + \lambda_2^2 + \lambda_3^2}}$$

- FA(D) is intrinsic to the tissue and is independent of fiber orientation, and diffusion sensitizing gradient directions
- FA(D) is useful to characterize the shape (degree of 'out-of-roundness') of the diffusion ellipsoid'
- Low FA:   High FA: 

Characterizing the Shape of the tensor: Fractional Anisotropy



Modules: Diffusion Tensor Scalar Measurements

3DSlicer

- Help & Acknowledgement
- Diffusion Tensor Scalar Measurements

Parameter set: Diffusion Tensor Scalar Measureme

Status

- IO

Input DTI Volume

Output Scalar Volume

- Operation

Estimation Parameters

- Trace
- Determinant
- RelativeAnisotropy
- FractionalAnisotropy
- Mode
- LinearMeasure
- PlanarMeasure
- SphericalMeasure
- MinEigenvalue
- MidEigenvalue
- MaxEigenvalue
- MaxFinenvalueProjectionX-

Default

• Data Probe

Red RAS: (69.3, 14.5, -0.8) Axial Sp: 1.5

L None()

F None()

Bdti (18, 54, 47) ColorOrientation 0

Set Input DTI Volume to 'dti'


Select Output Scalar Volume 'Create new Volume' and rename it 'fa'

Select the Operation 'Fractional Anisotropy'

Click on Apply to calculate the Fractional Anisotropy map of the tensor volume


Fractional Anisotropy

Select the Background volume to 'fa' in the red viewer and explore the FA values in the corpus callosum and in adjacent gray matter areas.



Fractional Anisotropy


Note how the FA values are high in the white matter areas, and low in gray matter regions




Fractional Anisotropy

Set the Foreground volume to 'None', and set the Background volume to 'dti' in the red viewer menu.


Go back to conventional layout





Part 2: Visualizing the tensor data

3D Visualization: Glyphs




3D Visualization: Glyphs

The screenshot shows the 3D Slicer software interface. The 'Modules' dropdown is set to 'Volumes'. The 'Active Volume' is 'dti'. The 'Display' panel is expanded, showing 'Scalar Display' with 'FractionalAnisotropy' selected. The 'Lookup Table' is 'Grey'. The 'Interpolate' checkbox is checked. The 'Window Level editor presents:' section shows 'W: 1.00', 'Auto W/L' (highlighted with a red arrow), and 'L: 0.50'. The 'Threshold' is set to 'Off'. The 'Glyphs on Slices Display' panel is expanded and highlighted with a red box, showing 'Slice Visibility' (Red, Yellow, Green), 'Opacity' (1.00), 'Scalar ColorMap' (Rainbow), 'Color by Scalar' (FractionalAnisotropy), and 'Scalar Range' (0.00 to 1.00). The 'Data Probe' is also visible. The main 3D view shows a brain slice with 'R', 'A', and 'L' labels. The bottom of the interface shows three orthogonal views (L, F, B) of the brain slice.

Click on Auto W/L to adjust the Window and Level values of the display

In the **Glyphs on Slices Display** panel, set the Color by Scalar parameter to 'ColorOrientation', and check Slice Visibility 'Red' '



3D Visualization: Glyphs



3D Visualization: Glyphs


Click on the link icon in the red slice viewer to unlink the three viewers.

Click on the eye icon to display the glyphs superimposed on the FA image in the 3D Viewer





L
F
B

3D Visualization: Glyphs



3D Visualization: Glyphs



Note the orientation of diffusion ellipsoid of the splenium of the corpus callosum (posterior part)



L
F
B

3D Visualization: Glyphs



Change the Glyph Type to 'Lines', and move the mouse inside the 3D viewer to refresh the display.



L
F
B

3D Visualization: Glyphs

Slicer displays the glyphs as lines that represent the principal direction of diffusion (main eigenvector)



L
F
B

3D Visualization: Glyphs

Select Red Slice Only layout in the layout menu

Window Level editor presets:

W: 1.00 Auto W/L L: 0.50

Threshold: Off
0.00 0.00

Histogram
Glyphs on Slices Display

Slice Visibility: ☒ Red ☐ Yellow ☐ Green

Opacity: 1.00

Scalar ColorMap: Rainbow

Color by Scalar: ColorOrientation

Scalar Range: 0 255

Glyph Type: Lines

Scale Factor: 18.00

Spacing: 15.00

Glyph EigenVector: Major

Data Probe

Red RAS: (-36.4, 27.1, -18.8) Axial Sp: 1.5

L None()

F dti (88, 46, 35) FractionalAnisotropy 0.07395

B dti (88, 46, 35) FractionalAnisotropy 0.07395

Change the Scale Factor and the Spacing and explore the glyphs in the optic chiasm area (slice S: -18.75)

Optic Chiasm

The optic chiasm corresponds to the part of the brain where the optic nerves cross.






Image from Gray's Anatomy

3D Visualization: Glyphs






Part 3: From tensors to tracts


Diffusion MRI tractography

- Tractography can be defined as the virtual reconstruction of the trajectory of water molecules along white matter bundles.
 - DTI tracts provide a mathematical representation of the underlying white matter anatomy.
 - Each voxel contains hundreds of thousands of axon fibers: size of a voxel $\sim 1\text{-}5\text{ mm}$ is very different from the diameter of an axon $\sim 0.1\text{-}10\text{ }\mu\text{m}$
- A DTI tract is not equivalent to a real fiber.



Tractography Seeding: ROI definition




ROI Definition



Set the Master Volume to 'fa'
Click on Apply in the pop-up window to create an empty labelmap 'fa-label'




ROI Drawing



Use the draw tool to outline the contour of the corpus callosum in the sagittal slice, and press Enter. Repeat the same operation on 3 adjacent sagittal slices.




ROI Drawing




Streamline tractography

Underlying Assumption: the orientation of the fibers is collinear with the direction of the principal eigenvector



Labelmap Seeding: I/O




Select the module **Tractography Label Map Seeding**

Set the Input DTI Volume to 'dti'
Set the Input Label Map to 'fa-label'

Set Output Fiber Bundle to 'Create New Fiber Bundle' and rename it 'corpusCallosum'

L
F
B

Labelmap Seeding: parameters



Select the Seed Placement Options to 'Use Index Space'.

Select Stopping Mode 'Fractional Anisotropy'

Select the default tractography Seeding parameters:

- Minimum length: 10 mm

- Maximum length: 800 mm


- Stopping value: 0.15

- Stopping track curvature: 0.8


- Integration step length: 0.5 mm

Click on **Apply**

Labelmap Seeding: Tracts





Labelmap Seeding: Tracts




Tractography Results

Select the Fiber Bundle corpusCallosum, and click on the eye icon next to Glyph




Slicer displays the glyphs (ellipsoids) along the tracts.


Tractography Results



Tractography Results




Tractography Results




Tractography Results



Tractography Results





Fiducial Seeding




Fiducial Seeding

Change the name of the fiducial to 'LeftCingulum'



Fiducial Seeding



Select the module
Tractography Interactive Seeding

Select the DTI volume 'dti'
Select the Fiducial List 'Fiducials List'
Select the Output Fiber Bundle 'Create New Fiber Bundle' and rename it 'Cingulum'

3D Slicer

Modules: Tractography Fiducial Seeding

Parameter set: FiducialSeedingParameters

IO

Input DTI Volume: Select a Di...ensorVolume

Input Fiducial List or Model: Fiducials List

Output Fiber Bundle: Select a FiberBundle

Seed Placement Options

Fiducial Region Size: 2.50mm

Fiducial Seeding Step Size: 1.00mm

Seed Selected Fiducials: ☐

Max Number of Seeds: 100

Tractography Seeding Parameters

Minimum Path Length: 20.00mm

Stopping Criteria: Fractional Anisotropy

Stopping Value: 0.25

Stopping Track Curvature: 0.70

Integration Step Length: 0.50mm

Enabling Options


Create Tracts Initially As: Tubes

Enable Seeding Tracts: ☒

Data Probe

L
F
B

Fiducial Seeding




Set the tractography parameters as follows:

- Fiducial region size: 2.5 mm
- Fiducial Seeding Step Size: 1.0 mm
- Maximum number of seeds: 100
- Minimum Path Length: 10 mm
- Stopping Criteria: Fractional Anisotropy
- Stopping Value: 0.15
- Stopping Track Curvature: 0.8
- Integration step length: 0.5 mm
- Create Tracts Initially as Tubes

Fiducial Seeding

Part of the left cingulum appears in the 3D viewer.

Move the Left Cingulum fiducial to explore the spatial relationship between the left cingulum and the corpus callosum



Fiducial Seeding

Click on the arrow icon to create a new fiducial, and position it in the right cingulum area.

Change the name of the new fiducial to 'Right Cingulum' in the Annotations module

3DSlicer

Modules: Annotations

Help & Acknowledgement

Annotations

Edit

Active list: Fiducials List

Vis	Lock	Edit	Value	Name	D
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		All Annot...	
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		Fiducials...	
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	-3.7, 12.0, 26.8	LeftCing...	
<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	8.5, 12.9, 24.9	RightCin...	

Data Probe

L
F
B

R S: 25.50

Y R: 12.00

G A: 13.50


RightCingulum

RightCingulum


Fiducial Seeding

Part of the left and right cingulum appear in the 3D viewer.


Move the fiducials to explore the spatial relationship between the left and right cingulum, and the corpus callosum




Fiducial Seeding



Fiducial Seeding




Tractography 'on-the-fly'




The Fiducial Seeding functionality allows you to do tractography 'on-the-fly' to explore white matter structures interactively

DTI Analysis



Conclusion



This tutorial guided you through the different steps of a Diffusion MR Analysis pipeline, from tensor estimation to 3D tracts visualization, for exploring and studying the brain white matter pathways.




Acknowledgments



National Alliance for Medical Image Computing
NIH U54EB005149



Neuroimage Analysis Center
NIH P41RR013218



Contact:
spujol@bwh.harvard.edu

Slicer Community

- www.slicer.org
- Mailing lists:
slicer-user@bwh.harvard.edu
slicer-devel@bwh.harvard.edu